
Trust Diagnostics



Overview

• Full DICOM Viewer + DICOM server, works as a PACS with single reporting workstation



Specifications

- DICOM server to archive DICOM studies from modalities
- Read and view DICOM files
- Read and view DICOM Directory
- Display images from different studies or examinations
- Displaying DICOM attributes of selected image
- · Playing cine DICOM with various speeds
- Query/Retrieve DICOM studies from PACS
- Brightness/contrast control
- Zooming and panning DICOM images
- · Reporting and voice notes for DICOM studies
- Export DICOM file to windows formats (jpg, etc)
- . Export cine DICOM file to windows formats (TIFF, AVI, etc)
- Make DICOM Directory of DICOM studies
- Print DICOM images via paper printers and DICOM printers
- Send DICOM studies to PACS